```
1 function spot results = FindTrajects(image label, start frame, end frame)
 2 %
 3 % Created by Isabel Llorente-Garcia, August 2011.
 4 % If you use this code please acknowledge Isabel Llorente-Garcia in your
 5 % publications.
 6 %
 7 % This function takes an image label such as '513', '490', etc...
 8 % which corresponds to a certain .sif image sequence in current folder,
 9 % and then finds all trajectories of the bright fluorescent spots in time.
10 %
11 % Example of how to use this function:
12 % for image "Heiko Thu Jun 24 2010 554.sif" in current folder:
13 % [numFrames2 frame Ysize2 frame Xsize2 image data2] = extract image sequence data('554');
14 %
15 % Reads .sif image sequence data and then for each frame:
16 % finds candidate spots,
17 % joins them to spots found and accepted on previous frame,
18 % eliminates coincidences in those candidates (points closer than 1 pixel),
19 % within all candidates, it finds actual spot centres through iterative Gaussian masking,
20 % and accepts only found spot centres with clipping flag equal to zero, constrained width and large enough SNR,
21 % eliminates coincidences from previously accepted found spot centres.
22 % Resulting final accepted spot centres are saved in the structure array spot final.
23 % Then link spots into trajectory segments: For second frame do differently
24 % and check only spots in 1st frame and compare to spots in second frame.
25 % For rest of frames (for frame k): A) first check loose spots (TrajNumber=0) two frames
26 % ago (k-2) and compare to spots in current frame (k).
27 % B) then check all spots in previous frame (k-1) and compare to spots in
28 % current frame (k).
29 % To decide on the best asignment of pairs of spots (link), we build up
30 % matrices of pair-wise distances, ratio of intensities and ratio of sigmas
31 % of all pairs of spots in the two frames being compared, and take the
32 % winning asignment as that with the smallest pairwise distance.
33 % The intensity ratio and ratio of sigmas has to be within certain min and
34 % max bounds too.
35 %
36 % start frame and end frames are the frames through which the loop runs to
37 % find centres of bright spots.
```

```
38 %
39 % Output:
40 % The output, spot results is a cell array with two elements:
41 % spot results = {params, spot final};
42 % The first element in the cell array, spot results{1}, contains all parameters used to run
43 % the function: "params" (this is a structure array itself).
44 % The second element in the cell array, spot results {2}, contains the track
45 % results: "spot final".
46 % "spot_final", is itself a structure array with end_frame by L elements,
47 % each of which is a structure with fields: CentreX, CentreY, IspTot,
48 % Sigma, IbqAvq, IbqTot, IinnerTot, ClipFlag, FrameNumber and SpotNumber.
49 % Along the dimension L, we have the different bright spots found on each
50 % frame (L will usually be the number of spot centres found in frame start,
51 % or the largest number of spots ever accepted).
52 % For example, spot final(100,8) is a structure with the above fields
53 % corresponding to the eighth found spot on frame 100.
54 %
55 % Example of params: spot_results{1}:
56 %
57 %
                 image label: '498'
                 start frame: 7
58 %
59 %
                   end frame: 500
60 %
                           r: 5
         max num candidates: 2000
61 %
62 %
          subarray halfwidth: 8
63 %
         inner circle radius: 5
64 %
            gauss mask sigma: 2
                sigmaFit min: 2
65 %
66 %
                sigmaFit max: 4
67 %
                     SNR min: 2
68 %
                     rsq min: 0.2000
69 %
                     deflate: 1
70 %
              d coincid cand: 1
             d coincid found: 1
71 %
72 %
                    d 01 max: 5
73 %
               Iratio 01 min: 0.5000
74 %
               Iratio 01 max: 3
```

```
75 %
            SigmaRatio 01 min: 0.5000
 76 %
            SigmaRatio 01 max: 2
 77 %
                    d 02 max: 5
 78 %
               Iratio 02 min: 0.5000
 79 %
               Iratio 02 max: 3
            SigmaRatio 02 min: 0.5000
 80 %
 81 %
           SigmaRatio 02 max: 2
 82 %
                         rei: 10000
 83 %
                   file name: 'Heiko Thu Jun 24 2010 498.sif'
 84 %
                   numFrames: 500
 85 %
                 frame Ysize: 512
 86 %
                 frame Xsize: 512
 87 %
 88 % Example of how to call this function:
 89 % s0 = FindTrajects('470',5,500); uses
 90 % image sequence 470 and finds trajectories in it, for frames 5 to 500.
 91 % Another example: s1 = FindTrajects('490',100,110);.
 92 % -----
 93 % Note: for reading .sif files you need 'read sif data direct.m',
 94 % 'GetAndorSifSize', etc. which are currently in path:
 95 % 'C:\Isabel\myMatlabFiles\IO Input\'.
 96 %
 97 % NOTE: before running this function you should move into the directory which
 98 % contains the image sequence data (labelled by 'image label').
 99
100
101 %% DEFINITIONS and PARAMETERS:
102
103 % I have saved all parameter values in the file "paramsForFindTrajects.m"
104 % in the current directory. Just calling the name of that file loads the
105 % parameter values into the workspace:
106 paramsForFindTrajects
107 % In this way, we can save different parameter sets for different data sets
108 % in an easy manner and run two Matlabs at the same time working with different parameter sets.
109
110 % %
111 % % define data (.sif image) directory:
```

```
112 % % dir data = 'Z:\Leake\Heiko Data\';
113 % % dir data = 'C:\Isabel\ExperimData\HeikoData\';
114 % % Print data directory on command window to guide user:
115 % disp(' ') % empty line
116 % disp(['The data directory (.sif images) is: ',cd]) % display current directory.
117 %
118 % % Save input parameters to "params":
119 % params.image label = image label;
120 % params.start frame = start frame;
121 % params.end frame = end frame;
122 %
123 % % Number of frames to average over when calculating a frame average in
124 % % order to get a Signal Mask to distinguish cell region from background
125 % % region:
126 % r = 5; % number of frames to average over, starting from start frame. (default: 5, end frame-start frame)
127 % % Save parameters to results as structure "params":
128 \% params.r = r;
129 %
130 % % If you don't want to use the cell mask to exclude spots found outside
131 % it, make the following parameter equal to 1 (then spots can be accepted anywhere on image):
132 % doNotUseCellMaskToTellSpotsRegion = 0;
133 % params.doNotUseCellMaskToTellSpotsRegion = doNotUseCellMaskToTellSpotsRegion;
134 %
135 % % Maximum number of candidate spots (if eq. 260000 candidate spots are
136 % % found, we get an error in function pdist: "Distance matrix has more
137 % % elements than the maximum allowed size in MATLAB"), hence, we limit the
138 % % max number of candidate spots:
139 % max num candidates = 2000; % should be around 200.
140 % % Save parameters to results as structure "params":
141 % params.max num candidates = max num candidates;
142 %
143 % PARAMETERS for finding spot centres (see findSpotCentrelframe.m, inputs to the function):
144 % % The total integrated spot intensity is a bqnd corrected one, inside a
145 % % circular mask of radius inner circle radius.
146 % subarray halfwidth = 8; % (Default: 8 pixels). Halfwidth of image square subarray
147 % % ROI, typically a square of size 17x17 pixels.
148 % inner circle radius = 5; % (Default: 5 pixels). Radius of inner circular mask that moves inside the fixed square subarray.
```

```
149 % gauss mask sigma = 2; % (Default: 2 pixels). Size in pixels of the applied Gaussian mask.
150 % quess sigma Fit = 3; % starting quess for Gaussian fit of brightspot intensity (default = 3).
151 % % Save parameters to results as structure "params":
152 % params.subarray halfwidth = subarray halfwidth;
153 % params.inner circle radius = inner circle radius;
154 % params.gauss mask sigma = gauss mask sigma;
155 % params.quess sigma Fit = quess sigma Fit;
156 %
157 % % PARAMETERS for deciding if we accept a spot centre found by the function
158 % % findSpotCentrelframe or not:
159 % sigmaFit min = -inner circle radius; % minimum acceptable sigma of gaussian fit to spot, in pixels (2) (-3).
160 % sigmaFit max = inner circle radius; % maximum acceptable sigma of gaussian fit to spot, in pixels (4) (3).
161 % SNR min = 1.4; % minimum acceptable signal-to-noise ratio (at least 2) (SNR as defined in findSpotCentrelframe.m).
162 % rsg min = 0.1; % minimum acceptable r-square value (0.2) (goodness of gaussian fit to spot).
163 % % Save parameters to results as structure "params":
164 % params.sigmaFit min = sigmaFit min;
165 % params.sigmaFit max = sigmaFit max;
166 % params.SNR min = SNR min;
167 % params.rsq min = rsq min;
168 %
169 % % PARAMETER to decide if deflation method is applied or not (subtract each
170 % % found spot to allow detection of weaker spots):
171 % deflate = 1;
172 % % Save parameters to results as structure "params":
173 % params.deflate = deflate;
174 %
175 % % PARAMETERS for eliminating coincident spots:
176 % d coincid cand = 1; % distance (in pixels) for eliminating coincidences in spot candidates.
177 % d coincid found = 1; % distance for eliminating coincidences in found spot centres.
178 % % Save parameters to results as structure "params":
179 % params.d coincid cand = d coincid cand; % distance for eliminating coincidences in spot candidates.
180 % params.d coincid found = d coincid found; % distance for eliminating coincidences in found spot centres.
181 %
182 % % PARAMETERS for building trajectories:
183 % % For linking spots in current and previous frames:
184 % d 01 max = 5; % max distance in pixels between spot centres in current and previous frames, for linking them into a ▶
trajectory (5).
```

```
185 % Iratio 01 min = 0.5; % min ratio of total spot intensities (after bond subtraction) (0.5).
186 % Iratio 01 max = 3; % max ratio of total spot intensities (after bgnd subtraction) (frame k-1/frame k) (large enough value ✓
(3) to account for blinking).
187 % SigmaRatio 01 min = 0.5; % min ratio of spot widths (sigma of Gaussian fit) (0.5).
188 % SigmaRatio 01 max = 2; % max ratio of spot width (sigma of Gaussian fit) (2).
189 % % Save parameters to results as structure "params":
190 % params.d 01 max = d 01 max;
191 % params. Iratio 01 min = Iratio 01 min;
192 % params. Iratio 01 max = Iratio 01 max;
193 % params. SigmaRatio 01 min = SigmaRatio 01 min;
194 % params. SigmaRatio 01 max = SigmaRatio 01 max;
195 %
196 % % For linking loose spots in current frame and 2 frames ago (jump of 1 frame in trajectory):
197 % d 02 max = 5; % max distance in pixels between spot centres in current frame and 2 frames ago. (default: 5)
198 % Iratio 02 min = 0.5; % min ratio of total spot intensities (after bgnd subtraction).
199 % Iratio 02 max = 3; % max ratio of total spot intensities (after bgnd subtraction).
200 % SigmaRatio 02 min = 0.5; % min ratio of spot widths (sigma of Gaussian fit).
201 % SigmaRatio 02 max = 2; % max ratio of spot width (sigma of Gaussian fit).
202 % % Save parameters to results as structure "params":
203 % params.d 02 max = d 02 max;
204 % params. Iratio 02 min = Iratio 02 min;
205 % params. Iratio 02 max = Iratio 02 max;
206 % params.SigmaRatio 02 min = SigmaRatio 02 min;
207 % params.SigmaRatio 02 max = SigmaRatio 02 max;
208 %
209 % % Use a very large number (larger than image size in pixels) for rejected asignments:
210 % rei = 10000;
211 % params.rej = rej; % Save parameters to results as structure "params".
212 %
213 % Parameter to exclude a region from accepting spots (see later on, lines 322 and 507):
214 % exclude_region = 0;
215 % % Note that when exclude region is 1, at the moment is it set for an image
216 % % with two channels, where spots are excluded also from around a horizontal
217 % % line centred on the image.
218 % exclude region width = subarray halfwidth; % exclude edges of image, only look at spots with a centre far enough from image &
edges.
219 % params.exclude region = exclude region;
```

```
220 % params.exclude region width = exclude region width;
221
222 % Alternative way of selecting image sequence file:
223 % uigetfile opens a file dialog box to choose data file:
224 % [file data,path data] = uigetfile({'*.sif'}, 'Chose image data sequence:');
225 % strcat('data (.sif image):',' ',path data,file data)
226 % open a file dialog box to choose analysis file:
227 % [file_analysis,path_analysis] = uigetfile({'*.xls'}, 'Chose analysis file (trajectory):');
228 % strcat('analysis file (.xls trajectory):',' ',path analysis,file analysis)
229
230 disp(' ') % empty line
231 disp(['The start frame for finding bright spot trajectories will be ',num2str(start frame)]) % start frame is an input.
232 disp(['The end frame for finding bright spot trajectories will be ',num2str(end frame)]) % end frame is an input.
234
235
236 %% Read in the image-sequence data:
237
238 % Read image-sequence file:
239 [numFrames frame Ysize frame Xsize image data image path] = extract image sequence data(image label);
240 % See "extract image sequence data.m".
241 % numFrames is the number of frames in the image sequence.
242 % To get frame number "p" do: image data(p).frame data.
243 % Frame dimensions are frame Ysize and frame Xsize.
245
246 % Save to parameters:
247 params.file name = image path;
248 params.numFrames = numFrames;
249 params.frame Ysize = frame Ysize;
250 params.frame Xsize = frame Xsize;
251
252
253 %% Calculate Signal Mask to distinguish cell region from background region:
254 % Use frame average (of first frames only) to calculate signal mask.
255
256 % Initialise frame accumulation in order to later calculate a frame average:
```

```
257 frame accumul = zeros(frame Ysize, frame Xsize);
258
259 % r is the number of frames to average over, starting from start frame.
260 % See PARAMETERS section.
261
262 for k = start frame:start frame+r % loop through frames.
263
        % Get frame data:
264
        frame = image data(k).frame data; % extract frame data, stored in the field 'frame data'.
265
        frame = double(frame);
        % Accummulate frames to then calculate frame average:
266
267
        frame accumul = frame accumul + frame;
268 end
269
270 % Calculate frame average as the accumulation of all frames divided by the number of frames:
271 frame avg = frame accumul/(r+1);
272
273 frame avq Gray = mat2gray(frame avq); % The input to function "getCellMaskAndBoundary" needs to be a grayscale image:
274
275 % Get SignalMask to know where cells are, to distinguish cells from background:
276 [SignalMask CellBoundaryMask] = getCellMaskAndBoundary(frame avg Gray);
277 % SignalMask is a matrix with 1 at positions where cells are and 0 at background.
278
279 % Or use getCellMaskAndBoundary2(frame,local region),
280 % with local region = [xleft xright ytop ybottom].
281 % [SignalMask CellBoundaryMask] = getCellMaskAndBoundary2(frame avg Gray,[1 frame Xsize round(frame Ysize/2) frame Ysize]); % 🗸
Use bottom half of image only.
282 % Using a local region for thresholding and finding the cell mask, makes this function much faster,
283 % since only spots within that cell mask will be considered.
284
285 % CHECK: uncomment the following line if you don't want to use a signal
286 % mask. Detected spots don't need to be within a signal mask:
287
288 if doNotUseCellMaskToTellSpotsRegion ==1
289
        SignalMask = ones(size(frame avg Grav.1), size(frame avg Grav.2));
290 end
291
292 %% Obtain candidate bright spots for start frame (first frame), and find spot centres for those:
```

```
293
294 frame = image data(start frame).frame data; % extract matrix data for first frame.
295 frame = double(frame);
296
297 disp(['frame number: ',num2str(start frame)]) % print frame number to Command Window.
298
299 % Xpos is a matrix of the same size as frame, containing x values for all
300 % pixels and similarly for Ypos (used in future sections):
301 [Xpos, Ypos] = meshgrid(1:frame Xsize, 1:frame Ysize);
302 % Note that the image thresholding occurrs in two halves: separating top and bottom halves.
303 % Find candidate-bright-spots on first frame:
304 frame Gray = mat2gray(frame); % The input to function "findCandidateSpots" needs to be a grayscale image:
305
306 [candidate spotsX 00 candidate spotsY 00] = findCandidateSpots(frame Gray,2); % Second input: use method 2, which seems to \( \mathbf{L} \)
work better.
307 % See C:\Isabel\myMatlabFiles\findCandidateSpots.m.
308 % candidate spotsX and candidate spotsY are two column vectors of the same
309 % length containing the x and y coordinates of the candidate bright spots found on the image.
310 % They contain integer numbers: coordinates or pixel numbers which give
311 % position on image plane.
312
313 % Reject candidate spots outside cell region (to speed up algorithm):
314 candidate spotsX 0 = []; % initialise empty vectors before loop.
315 candidate spotsY 0 = [];
316 for nn = 1:length(candidate spotsX 00)
       % Only use candidates for which there is a 1 in the SignalMask image:
317
318
       if SignalMask(candidate spotsY 00(nn), candidate spotsX 00(nn))==1
319
           candidate spotsX 0 = [candidate spotsX 0; candidate spotsX 00(nn)];
320
           candidate spotsY 0 = [candidate spotsY 0; candidate spotsY 00(nn)];
321
       end
322 end
323
324 disp(['no. of new candidate spots on start frame: ',num2str(length(candidate spotsX 0))])
325
326 % Error control:
327
        % Limit the max number of candidate spots (if eq. 260000 candidate spots are
328
        % found, we will get an error in function pdist: "Distance matrix has more
```

```
329
        % elements than the maximum allowed size in MATLAB").
330
        % Select only the first max num candidates then.
331
        if length(candidate spotsX 0) > max num candidates
332
            candidate spotsX 0 = candidate spotsX 0(1:max num candidates);
            candidate spotsY 0 = candidate spotsY 0(1:max num candidates);
333
            disp(['NOTE!! no. of candidate spots has been limited to ',num2str(max num candidates)])
334
335
        end
336
337 % % Check graphically:
338 % imshow(frame Gray,[]);
339 % hold on;
340 % plot(candidate spotsX 0, candidate spotsY 0, '*');
341 % figure;
342
343 % Find spot centres and decide if we accept them or not:
344 n =1; % Initialise index n (index for accepted spot centres which have a clipping flag equal to zero):
345 frame to search = frame; % Initialise frame to search for spot centres.
346
347 % Find spot centre through iterative masking:
348 for m = 1:size(candidate spotsX 0.1) % loop through all candidate spots.
        % Now find centre of bright spot using function findSpotCentrelframe:
349
350
        % use candidate spots as initial estimates and then iterate to find spot centre.
351
        % Image subarray ROI is a square of size 17x17 pixels (halfwidth is
        % 8 pixels), inner circular mask that moves inside the fixed 17x17
352
        % square has a radius of 5 pixels and the applied Gaussian
353
354
        % mask has a sigma of 2 pixels:
355
        spot result = findSpotCentrelframe(frame to search.candidate spotsX 0(m),candidate spotsY 0(m),subarray halfwidth, &
inner circle radius, gauss mask sigma, guess sigma Fit);
        spot result.FrameNumber = start frame; % Add new field containing frame number (time) to result structure.
356
357
358
359
        if (spot result.ClipFlag == 0 && spot result.noConverge == 0 && ...
360
                spot result.SigmaFit <= sigmaFit max && ...</pre>
                spot result.SigmaFit >= sigmaFit min && ...
361
362
                spot result.SNR >= SNR min &&...
                spot result.rsqFit >= rsq min) &&...
363
364
                (exclude region == 0 | | (exclude region ==1 && (spot result.CentreY <(frame Ysize/2-exclude region width) | | ✓
```

```
spot result.CentreY > (frame Ysize/2+exclude region width))))
365
                % Only accept and save result of found spot if clipping flag =0 and if values of sigmaFit, signal to noise and ✓
rsquare of fit are acceptable.
366
                spot result.SpotNumber = n; % Add new field containing spot number to result structure.
367
                spot final(start frame,n) = spot result; % store "good" found spots.
                % This is also saved in the final result spot final, structure array.
368
369
                % first index is for frame number, second index is for spot number.
370
                §_____
371
372
                if deflate==1 % see parameter section at the beginning.
373
                    % "Deflation" process: subtract from raw frame image the corresponding
374
                    % Gaussian fit of each found and accepted spot before finding next spot centre (enables acceptance of dimmer ⊌
spots).
375
376
                    % Matrices containing the x and y positions in the image frame: Xpos and Ypos
377
                    % Xpos is a matrix of the same size as frame, containing x values for all pixels and similarly for Ypos.
378
                    % Calculate Xpos, Ypos at the beginning: [Xpos, Ypos] = meshgrid(1:frame Xsize, 1:frame Ysize);
379
                    % Parameters of Gaussian fit of previously accepted spot:
                   x fit = spot final(start frame,n).CentreX;
380
                   v fit = spot final(start frame,n).CentreY;
381
382
                   I fit = spot final(start frame,n).IOFit;
383
                    sigma fit = spot final(start frame,n).SigmaFit;
                    % deflated frame (frame with found spot subtracted):
384
385
                    deflated frame = frame to search - I fit*exp(-((Xpos-x fit).^2+(Ypos-y fit).^2)/(2*sigma fit^2));
386
                    frame to search = deflated frame; % update frame to search for finding next spot-centre.
                    % % Graphical check of deflated frames:
387
388
                    % subplot(1,2,1); imshow(frame,[]); % frame is the original frame (always the same).
                    % subplot(1,2,2); imshow(deflated_frame,[]);
389
390
                end
391
392
393
                n = n+1; % advance index n for accepted spot centres.
394
        end
395 end
396
397 % % display the number of accepted spot-centres for this frame:
398 disp(['no. of accepted spot centres in first frame: ',num2str(n-1)])
```

```
399
400 % convert results of found spot-centre positions to a useful form that can
401 % be used as input candidate-spots on the following frame:
402 if (n-1) == 0 % error control: if no spots were accepted.
403
        found spot CentreX = [];
        found spot CentreY = [];
404
405
        % I need to create the whole spot final structure with all its fields
406
        % here, just in case the number of accepted spots in the first frame is
        % zero, in order not to get error: "Subscripted assignment between
407
        % dissimilar structures".
408
        % Save empty spot (we need this, otherwise if in the last frame the no. of accepted spots is 0, there will be no result ₺
409
spot final(end frame,:) and the following functions will fail).
        spot final(start frame,n).CentreX = [];
410
411
        spot final(start frame,n).CentreY = [];
412
        spot final(start frame,n).IspTot = [];
413
        spot final(start frame,n).rsgFit = [];
414
        spot_final(start_frame,n).SigmaFit = [];
415
        spot final(start frame,n).IOFit = [];
416
        spot final(start frame,n).bq noise offset afterBGsubtract = [];
417
        spot final(start frame,n).BqNoiseStd = [];
        spot final(start frame,n).IbqAvq = [];
418
419
        spot final(start frame,n).IbqTot = [];
420
        spot final(start frame,n).SNR = [];
        spot final(start frame,n).linnerTot = [];
421
422
        spot_final(start_frame,n).ClipFlag = [];
423
        spot_final(start_frame,n).noConverge = [];
424
        spot final(start frame,n).TraiNumber = [];
425
        spot final(start frame,n).FrameNumber = [];
426
        spot_final(start_frame,n).SpotNumber = [];
427 else
        found spot CentreX = [spot final(start frame,:).CentreX]'; % column vector with found CentreX positions of all candidate \( \sum \)
428
spots.
429
        found spot CentreY = [spot final(start frame,:).CentreY]'; % column vector with found CentreY positions of all candidate \( \mathbb{E} \)
spots.
430 end
431
432 % Check graphically:
```

```
433 figure;
434 imshow(frame Gray,[]);
435 hold on;
436 plot(found spot CentreX, found spot CentreY, 'o', 'Color', 'q', 'MarkerSize', 10) % plot accepted spot centres in green.
437 pause(0.1); % this pause is needed to give time for the plot to appear
438 hold off;
440
441
442 %% Loop through selected frames:
443
444 tr =1; % initialise trajectory index.
445
446 for k = (start frame+1):end frame
447
        % to go through all frames do instead: for k = 1:length(sifData)
448
449
        frame = image data(k).frame data; % extract frame data which is stored in field 'frame data'.
450
        frame = double(frame);
451
452
        imshow(frame,[],'Border','tight','InitialMagnification',150); % show image scaled between its min and max values ([]).
       hold on;
453
454
455
        disp(['frame number: ',num2str(k)]) % print frame number to Command Window.
456
        &_____
457
        % Find new candidate spots for this frame:
458
459
        frame Gray = mat2gray(frame); % The input to function "findCandidateSpots" needs to be a grayscale image:
460
461
        [candidate spotsX 00 candidate spotsY 00] = findCandidateSpots(frame Gray,2); % Second input: use method 2, which seems 
to work better.
462
        % the subindex " 00" in candidate spotsX 00 indicates newly found spot
463
        % candidates for the current frame. On the other hand, found spot CentreX and found spot CentreY are the
464
        % accepted spot-centre positions coming from the previous frame.
465
466
        % Reject candidate spots outside cell region (SignalMask) (to speed up algorithm):
467
        candidate spotsX 0 = []; % initialise empty vectors before loop.
468
        candidate spotsY 0 = [];
```

```
469
        for nn = 1:length(candidate spotsX 00)
470
            % Only use candidates for which there is a 1 in the SignalMask image:
471
            if SignalMask(candidate spotsY 00(nn),candidate spotsX 00(nn))==1
472
                candidate spotsX 0 = [candidate spotsX 0; candidate spotsX 00(nn)];
473
                candidate spotsY 0 = [candidate spotsY 0; candidate spotsY 00(nn)];
474
            end
475
        end
476
        disp(['no. of new candidate spots: ',num2str(length(candidate spotsX 0))])
477
478
479
480
        % Join accepted spot-centre positions from previous frame with
481
        % candidate spots for this frame to use them as new candidates for this frame:
482
        candidate spotsX = [found spot CentreX; candidate spotsX 0];
483
        candidate spotsY = [found spot CentreY; candidate spotsY 0];
484
485 %
               % Plot new candidate spots in vellow and found spot centres from
486 %
               % previous frame in cyan
                plot(candidate spotsX 0,candidate spotsY 0,'+','Color','y','MarkerSize',3);
487 %
488 %
                pause(0.5);
489 %
                plot(found spot CentreX, found spot CentreY, '+', 'Color', 'c', 'MarkerSize', 3);
490 %
                pause(0.5);
491
        disp(['no. of initial total candidate spots: '.num2str(length(candidate spotsX))])
492
493
494
        % Error control:
495
        % Limit the max number of candidate spots (if eq. 260000 candidate spots are
        % found, we get an error in function pdist: "Distance matrix has more
496
497
        % elements than the maximum allowed size in MATLAB").
498
        % Select only the first max num candidates then.
499
        if length(candidate spotsX) > max num candidates
500
            candidate spotsX = candidate spotsX(1:max num candidates);
501
            candidate spotsY = candidate spotsY(1:max num candidates);
            disp(['NOTE!! no. of initial total candidate spots has been limited to '.num2str(max num candidates)])
502
503
        end
504
505
```

```
506
        % Eliminate coincidences in spot candidates:
507
        [candidate spotsX candidate spotsY pos] = eliminateCoincidentSpots(candidate spotsX, candidate spotsY, d coincid cand);
508
509
        % see C:\Isabel\myMatlabFiles\eliminateCoincidentSpots.m
510
        % The function checks the distances between all pairs of points with x
        % and y coordinates candidate spotsX and candidate spotsY respectively,
511
512
        % and removes those points (x,y) which are closer than one pixel (distance<1) to
513
        % another point in the list.
514
515
              % for debugging:
516
              x2 = candidate spotsX;
517
        ે
              y2 = candidate spotsY;
518
        disp(['no. of total candidate spots after eliminating coincidences: ',num2str(length(candidate spotsX))])
519
520
521
              % Plot all candidate spots in magenta after removing coincidences:
522
               plot(candidate spotsX,candidate spotsY,'+','Color','m','MarkerSize',3);
523
               pause(0.5);
524
525
526
527
        n =1; % Initialise index n (index for accepted spot centres which have a clipping flag equal to zero):
528
        frame to search = frame; % Initialise frame to search for spot centres.
529
        for m = 1:size(candidate spotsX,1) % for each frame, loop throuh all the candidate spots.
530
531
            % Now find centre of bright spot using function findSpotCentrelframe:
            % use candidate spots as initial estimates and then iterate to find spot centre.
532
            % Image subarray ROI is a square of size 17x17 pixels (halfwidth is
533
            % 8 pixels), inner circular mask that moves inside the fixed 17x17
534
            % square has a radius of 5 pixels and the applied Gaussian
535
            % mask has a sigma of 2 pixels:
536
537
            spot result = findSpotCentre1frame(frame to search, candidate spotsX(m), candidate spotsY(m), subarray halfwidth, &
inner circle radius, gauss mask sigma, guess sigma Fit);
            % index k is for frame number, index m is for spot number
538
            spot result.FrameNumber = k; % Add new field containing frame number (time) to result structure.
539
540
541
            % accepted spot centres:
```

```
542
            if (spot result.ClipFlag == 0 && spot result.noConverge == 0 && ...
543
                    spot result.SigmaFit <= sigmaFit max && ...</pre>
                    spot result.SigmaFit >= sigmaFit min && ...
544
                    spot result.SNR >= SNR min &&...
545
546
                    spot result.rsqFit >= rsq min)&&...
547
                    (exclude region == 0 | (exclude region ==1 && (spot result.CentreY <(frame Ysize/2-exclude region width) | | \(\mu\)
spot result.CentreY >(frame Ysize/2+exclude region width))))
548
               % Only accept and save result of found spot if clipping flag =0 and if values of sigmaFit, signal to noise and ✓
rsquare of fit are acceptable.
               spot(k,n) = spot result; % store accepted found spots in this preliminary result.
549
550
               % first index is for frame number, second index is for spot number.
551
552
                           plot(spot(k,n).CentreX,spot(k,n).CentreY,'o','Color','r','MarkerSize',10); % Plot found centre spots ✓
in red
                           553
554
               if deflate==1 % see parameter section at the beginning.
555
                   % "Deflation" process: subtract from raw frame image the corresponding
                    % Gaussian fit of each found and accepted spot before finding
556
                    % next spot centre (enables acceptance of dimmer spots).
557
                   % Xpos and Ypos are matrices of x and y positions on image frame.
558
                    % Parameters of Gaussian fit of previously accepted spot:
559
560
                   x_{fit} = spot(k,n).CentreX;
561
                   v fit = spot(k,n).CentreY;
                   I fit = spot(k,n).IOFit;
562
                    sigma_fit = spot(k,n).SigmaFit;
563
                   % deflated frame (frame with found spot subtracted):
564
565
                   deflated frame = frame to search - I fit*exp(-((Xpos-x fit).^2+(Ypos-y fit).^2)/(2*sigma fit^2));
                    frame to search = deflated frame; % update frame for finding next spot-centre.
566
                   % % Graphical check of deflated frames:
567
                    % subplot(1,2,1); imshow(frame,[]); % frame is the original frame (always the same).
568
                    % subplot(1,2,2); imshow(deflated_frame,[]);
569
570
               end
571
572
573
               n = n+1; % and advance index n for accepted spot centres.
574
            end
575
```

```
576
        end
577
578
        % display the number of accepted spot-centres for each frame:
579
        disp(['no. of accepted spot centres: '.num2str(n-1)])
580
581
        % The following two lines are used together with the previous two
582
        % "plot" and "imshow" (commented off) lines:
583
        pause(0.1); % this pause is needed to give time for the plot to appear
584
            hold off;
585
        % convert results of found spot-centre positions to a useful form that can
586
587
        % be used as input candidate-spots on the following frame:
        if (n-1) == 0 % error control: if no spots were accepted.
588
589
            found spot CentreX = [];
590
            found spot CentreY = [];
591
            spot final(k,n). SpotNumber = []; % Save empty spot (we need this, otherwise if in the last frame the no. of accepted \checkmark
spots is 0, there will be no result spot final(end frame,:) and the following functions will fail).
592
        else
593
            found spot CentreX = [spot(k,:).CentreX]'; % column vector with found CentreX positions of all candidate spots.
594
            found spot CentreY = [spot(k,:).CentreY]'; % column vector with found CentreY positions of all candidate spots.
595
596
597
            % Eliminate coincidences in result of last found spots for a given frame (for distance <1):
598
            [found spot CentreX found spot CentreY pos final] = eliminateCoincidentSpots(found spot CentreX, found spot CentreY, &
d coincid found);
599
            % see C:\Isabel\myMatlabFiles\eliminateCoincidentSpots.m
600
            % pos final contains positions of selected, kept spot centres.
601
602
            % Save final spots to variable final spots:
603
            n=1; % index for final kept spot.
604
605
            for ii = 1:length(pos_final)
                mientras = spot(k,pos final(ii)); % intermediate result.
606
                mientras.SpotNumber = n; % Add new field containing spot number to result structure.
607
                spot final(k,n)=mientras; % final result structure of accepted spot centres.
608
                n = n+1;
609
610
            end
```

```
611
            % Plot found spot centres:
612
            pause(0.5);
613
            plot(found spot CentreX, found spot CentreY, 'o', 'Color', 'g', 'MarkerSize', 10) % plot final accepted spot centres in ✓
green.
614
            pause(0.1); % this pause is needed to give time for the plot to appear
            hold off;
615
616
617
            disp(['no. of final found spot centres after eliminating coincidences: '.num2str(length(found spot CentreX))])
618
        end
619
620
621
622
        %______
623
        % LINKING SPOTS INTO TRAJECTORY SEGMENTS:
624
625
        % Link found and accepted spots into trajectory segments:
626
627
        % Trajectory index tr is initialised to 1 outside the loop through frames (k loop).
628
629
        % Do differently FOR SECOND FRAME (k == start frame+1): compare only accepted spots in
        % previous and current frames:
630
631
        if k == start_frame+1 && ... % If second frame and
632
                (n-1)~=0 && ... % if the number of accepted spot centres is not zero and
                isempty([spot final(k-1,:).SpotNumber]) == 0 && ... % at least 1 accepted spot in previous frame and
633
                isempty([spot final(k,:).SpotNumber]) == 0 % at least 1 accepted spot in current frame.
634
635
            % There are no trajectories jet, so compare accepted spots in previous and current frames:
636
            NO = max(cat(1,spot final(k-1,:).SpotNumber)); % no. of accepted spots in previous frame.
637
            % Note: cat(1,spot final(k-1,:).SpotNumber) gives a column vector with the values of SpotNumber for all non-empty ✓
accepted spots in frame k-1.
            N1 = \max(\text{cat}(1, \text{spot final}(k, :).\text{SpotNumber})); % no. of accepted spots in current frame.
638
639
640
            % Create cell arrays with empty elements to pre-asign sizes:
641
            d01 = cell(N0,N1); % Note: d01 is a cell array (matrix) but d 01 below is a scalar.
            Iratio01 = cell(N0.N1); % Note: Iratio01 is a cell array (matrix) but Iratio 01 below is a scalar.
642
            SigmaRatio01 = cell(N0,N1); % Note: SigmaRatio01 is a cell array (matrix) but SigmaRatio 01 below is a scalar.
643
644
645
            for g0 = 1:N0 % loop though accepted spots in previous frame.
```

```
646
                for g1 = 1:N1 % loop though accepted spots in current frame.
647
                    % d 01: distance between spot centres in previous and current frames:
648
                    d 01 = sgrt((spot final(k-1,q0).CentreX-spot final(k,q1).CentreX)^2+(spot final(k-1,q0).CentreY-spot final(k, <math>\checkmark
q1).CentreY)^2);
649
                    % Iratio 01: ratio of intensities of spot centre in previous and current frames:
                    Iratio 01 = spot final(k-1,q0).IspTot/spot final(k,q1).IspTot;
650
651
                    % SigmaRatio 01: ratio of widths of spots (Gaussian fits) in previous and current frames:
652
                    SigmaRatio 01 = spot final(k-1,q0).SigmaFit/spot final(k,q1).SigmaFit;
653
654
                                       d 01
655
                                       Iratio 01
656
                                       SigmaRatio 01
657
658
                    % Accept and save trajectory if spots in previous and
659
                    % current frames fulfill the following conditions:
660
                    if d 01 < d 01 max && ... % see PARAMETERS at start of this function.
661
                            Iratio 01 min <= Iratio 01 && Iratio 01 <= Iratio 01 max && ...</pre>
                            SigmaRatio 01 min <= SigmaRatio 01 && SigmaRatio 01 <= SigmaRatio 01 max
662
                        % Asign accepted values to cell array elements to store them:
663
664
                        d01\{q0,q1\} = d 01; % use {} for cell arrays.
                        Iratio01{q0,q1} = Iratio 01;
665
666
                        SigmaRatio01{q0,q1} = SigmaRatio_01;
667
                    else % rejected asignments:
                        d01{q0,q1} = rej; % Use rej for asignments not accepted (images usually 512x512arrays, so rej pix is an ⊌
668
impossibly large distance, this is why it is chosen here).
                        Iratio01\{q0,q1\} = rej; % Use rej for asignments not accepted.
669
670
                        SigmaRatio01{q0,q1} = rej; % Use rej for rejected asignments.
671
                    end
672
                end
673
674 %
                              d01
                              [d01{q0,:}]
675 %
                % Note that [d01{q0,:}] gives only non-empty elements of row q0
676
                % in the cell array d01 as a row vector, that's why we had to
677
                % give a numeric value rej to non-accepted asignments.
678
679
680
                % Note that if all asignments in previous step are rejected.
```

```
% [d01{q0.:}]  will be a list of rej values, and its minimum will
681
682
                % be rej.
                % If list of "linkable" spots, [d01{q0,:}], has no accepted
683
                % asignments (all values are rei):
684
                if min([d01{q0,:}]) == rej
685
686
                    % Asign trajectory number 0 to the spot in the previous frame only:
687
                    spot final(k-1,q0).TrajNumber = 0;
688
                else % if there is at least one accepted asignment for a given spot in the previous frame:
689
690
                    % Decide of all possible accepted spots (in current frame)
691
                    % that could be linked to spot q0 in previous frame, which one is the best:
692
                    % We take the best as the closest one to spot q0:
                    gl chosen = find([d01{q0},:]) == min([d01{q0},:])); % find position of the minimum pair-wise distance.
693
694
695 %
                                   q1 chosen
696
697
                    % Check if there is a better competing asignment for a given spot gl in the current
                    % frame from another spot in the previous frame.
698
                    % Hence, check also column-wise in matrix d01 to avoid asigning a traj
699
700
                    % number to a spot of in the current frame that had already
701
                    % had a trai number asigned to it linking it to a different spot g0 in
702
                    % the previous frame, which might be at a shorter distance
703
                    % from it than the current one.
704
                                   [d01{:,ql chosen}] % chosen column of d01 matrix of distances.
705 %
706
707
                    % Asign trajectory numbers to structure spot final:
                    % If the found distance in that column is not the minimum one:
708
                    if q0 \sim find([d01{:,q1 chosen}] = min([d01{:,q1 chosen}]));
709
                        spot final(k-1,q0). TrajNumber = 0; % asign trajectory number 0 to spot in previous frame.
710
711
                    else
712
                        spot final(k-1,q0). TrajNumber = tr; % asign trajectory number to spot in previous frame, to spot final \checkmark
structure.
713
                        spot final(k,ql chosen). TraiNumber = tr; % asign same trajectory number to spot in current frame.
                        tr = tr+1; % advance trajectory-number index.
714
715
                    end
716
                end
```

```
717
            end
718
719
720
721
        else % for FRAMES k \ge start frame+2, from third chosen frame on:
722
723
            % A) Compare loose spots (TrajNumber is 0) two frames ago (k-2)
724
            % to found spots in current frame (TraiNumber is []):
725
            726
            % DO maybe!!
727
            if k \sim start frame+1 & (n-1)\sim 0 & \dots % If the number of accepted spot centres is not zero and
728
                    isempty([spot final(k-2,:).SpotNumber]) == 0 && ... % at least 1 accepted spot 2 frames ago and
729
                    isempty([spot final(k,:).SpotNumber]) == 0 % at least 1 accepted spot in current frame.
730
731
                N0 = \max(\text{cat}(1, \text{spot final}(k-2,:).\text{SpotNumber})); % no. of accepted spots 2 frames ago.
732
                N1 = \max(\text{cat}(1, \text{spot final}(k, :).\text{SpotNumber})); % no. of accepted spots in current frame.
733
734
                % Create cell arrays with empty elements to pre-asign sizes.
735
                d02 = cell(N0,N1); % Note: d02 is a cell array (matrix) but d 02 below is a scalar.
736
                Iratio02 = cell(N0.N1); % Note: Iratio02 is a cell array (matrix) but Iratio 02 below is a scalar.
                SigmaRatio02 = cell(N0.N1); % Note: SigmaRatio02 is a cell array (matrix) but SigmaRatio 02 below is a scalar.
737
738
739
                for g0 = 1:N0 % loop though loose accepted spots 2 frames ago.
740
                    if spot final(k-2,g0). TraiNumber == 0 % only for loose (unlinked) spots (TraiNumber=0) two frames ago (so \( \sigma \)
only rows q0 in matrix d01 which have unlinked spots will fill up).
741
                        for q1 = 1:N1 % loop though accepted spots in current frame.
742
                            % d 01: distance between spot centres in previous and current frames:
743
                            d 02 = sqrt((spot final(k-2,q0).CentreX-spot final(k,q1).CentreX)^2+(spot final(k-2,q0).CentreY-<math>\checkmark
spot_final(k,q1).CentreY)^2);
744
                             % Iratio 01: ratio of intensities of spot centre in previous and current frames:
745
                            Iratio 02 = spot final(k-2,q0).IspTot/spot final(k,q1).IspTot;
746
                            % SigmaRatio 01: ratio of widths of spots (Gaussian fits) in previous and current frames:
747
                            SigmaRatio 02 = \text{spot final}(k-2,q0).\text{SigmaFit/spot final}(k,q1).\text{SigmaFit};
748
749
                            9
                                                       d 02
750
                            응
                                                       Iratio 02
751
                                                       SigmaRatio 02
```

```
752
753
                            % Accept and save trajectory if spots in previous and
754
                            % current frames fulfill the following conditions:
755
                            if d 02 < d 02 max && ... % see PARAMETERS at start of this function.
756
                                     Iratio 02 min <= Iratio 02 && Iratio 02 <= Iratio 02 max && ...
                                     SigmaRatio 02 min <= SigmaRatio 02 && SigmaRatio 02 <= SigmaRatio 02 max
757
758
                                 % Asign accepted values to cell array elements to store them:
759
                                d02\{q0,q1\} = d 02; % use {} for cell arrays.
                                Iratio02\{q0,q1\} = Iratio 02;
760
761
                                SigmaRatio02{q0,q1} = SigmaRatio 02;
762
                            else % rejected asignments:
763
                                d02{q0,q1} = rej; % Use rej for asignments not accepted (images usually 512x512arrays, so rej pix ✓
is an impossibly large distance, this is why it is chosen here).
764
                                Iratio02\{q0,q1\} = rej; % Use rej for asignments not accepted.
                                SigmaRatio02{q0,q1} = rej; % Use rej for rejected asignments.
765
766
                            end
767
                        end
768
                                               d02
769 %
770 %
                                               [d02{q0,:}]
771
                        % Note that [d02\{q0,:\}] gives only non-empty elements of the cell array d01 as a row vector.
772
773
                        % Note that if all asignments in previous step are rejected,
                        % [d02{q0,:}] will be a list of rej values, and its minimum will be rej.
774
775
                        % If list of "linkable" spots, [d02\{q0,:\}], has no accepted asignments (all values are rej):
776
                        if min([d02{q0,:}]) == rej
777
                            % Asign trajectory number 0 to the spot two frames ago only:
                            spot final(k-2,q0).TrajNumber = 0; % point stays loose (unlinked).
778
779
780
                        else % if there is at least one accepted asignment for a given spot two frames ago:
781
782
                            % Decide of all possible accepted/saved spots (in current frame)
783
                            % that could be linked to spot q0 in previous frame (of all possible asignments), which one is the
best:
784
                            % We take the best as the closest one to spot q0:
785
                            q1 chosen = find([d02{q0},:]) == min([d02{q0},:])); % find position of the minimum pair-wise distance.
786
```

```
787 %
                                                    al chosen
788
789
                             % Check if there is a better competing asignment for a given spot gl in the current
790
                             % frame from another spot q0 two frames ago.
791
                             % Hence, check also column-wise in d01 to avoid asigning a traj
792
                             % number to a spot gl in the current frame that had already
793
                             % had a traj number asigned to it linking it to a different spot g0 two frames ago which might be at \checkmark
a shorter distance
794
                             % from it than the current one.
795
796 %
                                                     [d02{:,q1 chosen}] % chosen column of d01 matrix of distances.
797
798
                             % Asign trajectory numbers to structure spot final:
799
                             % If the found distance in that column is not the minimum one:
                             if q0 \sim find([d02\{:,q1 \text{ chosen}\}] = min([d02\{:,q1 \text{ chosen}\}]));
800
                                 % Asign trajectory number 0 to the spot two frames ago only:
801
802
                                 spot final(k-2,q0).TrajNumber = 0; % point stays loose (unlinked).
803
                             else
804
                                  spot final(k-2,q0). TrajNumber = tr; % asign trajectory number to spot two frames ago, to ✓
spot final structure.
805
                                  spot final(k,ql chosen). TraiNumber = tr; % asign same trajectory number to spot in current frame.
806
                                 tr = tr+1; % advance trajectory-number index.
807
                             end
808
                         end
809
                     end
810
                end
811
            end
812
813
            % B) Compare loose spots (TraiNumber is []) and trajectories (TraiNumber is >0)
814
815
            % in previous frame (k-1) to found spots in current frame:
816
            if (n-1)\sim=0 && ... % If the number of accepted spot centres is not zero and
817
                     isempty([spot final(k-1,:).SpotNumber]) == 0 \&\& ... % at least 1 accepted spot in previous frame.
                     isempty([spot final(k,:).SpotNumber]) == 0 % at least 1 accepted spot in current frame.
818
819
820
                % ZZZZZZZZZZZZZZZZZZZZZZZZ
821
                N0 = \max(\text{cat}(1, \text{spot final}(k-1,:).\text{SpotNumber})); % no. of accepted spots in previous frame.
```

```
822
                N1 = max(cat(1,spot_final(k,:).SpotNumber)); % no. of accepted spots in current frame.
823
824
                % Create cell arrays with empty elements to pre-asign sizes.
825
                d01 = cell(N0.N1);
                Iratio01 = cell(N0,N1);
826
                SigmaRatio01 = cell(N0,N1);
827
828
829
                for g0 = 1:N0 % loop though accepted spots in previous frame.
                    for g1 = 1:N1 % loop though accepted spots in current frame.
830
831
                        % d 01: distance between spot centres in previous and current frames:
832
                        d 01 = sgrt((spot final(k-1,q0).CentreY-spot final(k,q1).CentreX)^2+(spot final(k-1,q0).CentreY-&
spot_final(k,q1).CentreY)^2);
833
                        % Iratio 01: ratio of intensities of spot centre in previous and current frames:
                        Iratio 01 = spot final(k-1,q0).IspTot/spot final(k,q1).IspTot;
834
835
                        % SigmaRatio 01: ratio of widths of spots (Gaussian fits) in previous and current frames:
836
                        SigmaRatio 01 = spot final(k-1,q0).SigmaFit/spot final(k,q1).SigmaFit;
837
838
                        응
                                               d 01
839
                        ્ટ
                                               Iratio 01
840
                        응
                                               SigmaRatio 01
841
842
                        % Accept and save trajectory if spots in previous and
843
                        % current frames fulfill the following conditions:
                        if d 01 < d 01 max && ... % see PARAMETERS at start of this function.
844
                                Iratio 01 min <= Iratio 01 && Iratio 01 <= Iratio 01 max && ...</pre>
845
                                SigmaRatio 01 min <= SigmaRatio 01 && SigmaRatio 01 <= SigmaRatio 01 max
846
847
                            % Asign accepted values to cell array elements to store them:
                            d01\{q0,q1\} = d_01; % use {} for cell arrays.
848
                            Iratio01\{q0,q1\} = Iratio_01;
849
850
                            SigmaRatio01{q0,q1} = SigmaRatio 01;
851
                        else % rejected asignments:
852
                            d01{q0,q1} = rej; % Use rej for asignments not accepted (images usually 512x512arrays, so rej pix is ✓
an impossibly large distance, this is why it is chosen here).
                            Iratio01{q0,q1} = rej; % Use rej for asignments not accepted.
853
                            SigmaRatio01{q0,q1} = rej; % Use rej for rejected asignments.
854
855
                        end
856
                    end
```

```
857
858 %
                                      d01
                                      [d01\{q0,:\}] % last row of d01 matrix of distances.
859 %
                    % Note that [d01\{q0,:\}] gives only non-empty elements
860
                    % of row q0 in the cell array d01, as a row vector.
861
862
863
                    % Note that if all asignments in previous step are rejected,
864
                    % [d01{q0.:}] will be a list of rej values, and its minimum will
865
                    % be rej.
                    % If list of "linkable" spots, [d01{q0,:}], has no accepted
866
867
                    % asignments (all values are rej):
868
                    if min([d01{q0,:}]) == rej
869
870
                        if isempty(spot final(k-1,q0).TrajNumber) % if point in previous frame was not part of a trajectory ✓
(TrajNumber=[]):
871
                            % Asign trajectory number 0 to the spot in the previous frame only:
872
                            spot_final(k-1,q0).TrajNumber = 0;
873
                        end
874
875
                    else % if there is at least one accepted asignment for a given spot two frames ago:
876
877
                        % Decide of all possible accepted/saved spots (in current frame)
878
                        % that could be linked to spot q0 in previous frame, which one is the best:
879
                        % We take the best as the closest one to spot q0:
                        q1 chosen = find([d01{q0,:}]) == min([d01{q0,:}])); % find position of the minimum pair-wise distance.
880
881
882 %
                                          al chosen
883
                        % Check if there is a better competing asignment for a given spot g1 in the current
884
                        % frame from another spot q0 in the previous frame.
885
886
                        % Hence, check also column-wise in d01 to avoid asigning a traj
887
                        % number to a spot ql in the current frame that had already
                        % had a traj number asigned to it linking it to a different spot g0 in
888
                        % the previous frame which might be at a shorter distance
889
                        % from it than the current one.
890
891
892 %
                                          [d01{:,gl chosen}] % chosen column of d01 matrix of distances.
```

```
893
894
                        % If the found distance in that column is not the minimum one:
895
                        if q0 \sim find([d01{:,q1 chosen}] = min([d01{:,q1 chosen}]));
896
                            if isempty(spot final(k-1,q0).TrajNumber) % if point in previous frame was not part of a trajectory ✓
(TrajNumber=[]):
897
                                % Asign trajectory number 0 to the spot in the previous frame only:
898
                                spot final(k-1,q0).TrajNumber = 0;
899
                            end
900
                        else
901
                            % Asign trajectory numbers to structure spot final:
902
                            if spot final(k-1,q0). TrajNumber > 0 % if point in previous frame was already part of a trajectory:
903
                                 spot final(k,q1 chosen). TrajNumber = spot final(k-1,q0). TrajNumber; % asign that trajectory ✓
number to spot in current frame.
904
                            else % if point in previous frame was not part of a trajectory:
905
                                 spot final(k-1,q0). TrajNumber = tr; % asign new trajectory number to spot in previous frame.
                                 spot final(k,q1 chosen). TrajNumber = tr; % asign same trajectory number to spot in current frame.
906
907
                                tr = tr+1; % advance trajectory-number index.
908
                            end
909
                        end
910
                    end
911
                end
912
                % ZZZZZZZZZZZZZZZZZZZZZZZZ
913
            end
914
        end
915
916
917
918 end % loop through selected frames
919
920
921
922 %% OUTPUT OF SPOT-FINDING PROCESS: final output spot results:
923 %
924 spot results = {params, spot final};
925 % params is a structure array containing all parameters used to run the
926 % function.
927 % spot final, is a structure array with end frame x L elements.
```

```
928 % each of which is a structure with fields:
929 %
          'CentreX'
930 %
          'CentreY'
931 %
          'IspTot'
932 %
          'rsqFit'
933 %
          'SiqmaFit'
934 %
          'IOFit'
935 %
          'BaNoiseStd'
936 %
          'IbqAvq'
937 %
          'IbqTot'
938 %
          'SNR'
939 %
          'IinnerTot'
940 %
          'ClipFlag'
941 %
          'TrajNumber'
942 %
          'FrameNumber'
943 %
          'SpotNumber'
944 %
945 % Along the dimension L, we have the different bright spots found on each
946 % frame (L will often be the number of spot centres found in frame start,
947 % it is always the largest number of spots ever accepted on one frame).
948 % For example, spot final(100.8) is a structure with the above fields
949 % corresponding to the eighth found spot on frame 100.
950 %
951 % Note that even if we only analyse from start frame to end frame.
952 % spot final is a list containing empty structure arrays from index 1 to
953 % index start frame, and then the found spots for the analysed frames start frame to end frame.
954 %
955 % The result is padded to a fixed number of spot structures for reach
956 % frame (the maximum no. of accepted found spots of all frames), so that for a given
957 % frame in which less found spots have been accepted, the remaining
958 % elements are padded with empty structures with empty fields [].
959 % To check if a given spot is empty: isempty(spot final(101,10).CentreX)
960 % gives 1 if field "CentreX" of tenth spot found and accepted in frame 101
961 % is empty (equal to []).
962 %
963 % e.g. cat(1,spot final(100,:).SigmaFit) gives a vector column with all the
964 % non-empty SigmaFit values for all spot centres accepted in frame 100.
```